



东莞中力仪器科技有限公司  
Zhongli Instrument Technology Co., Ltd

---

**ZL-5060A**

**Protective clothing synthetic blood penetration tester**



Tel: [+86-769-3365 7986](tel:+86-769-33657986)  
Website: [www.zltester.com](http://www.zltester.com)

Fax: [+86-769-3365 7987](tel:+86-769-33657987)  
E-mail: [zhongli@zltester.com](mailto:zhongli@zltester.com)



# 东莞中力仪器科技有限公司


## Zhongli Instrument Technology Co., Ltd

**Apply for :** Used to test the resistance of synthetic medical clothing made of fabrics, nonwovens and other materials to the penetration of synthetic blood under different levels of test pressure.

**Standard:** GB19082; ISO16603; ASTM F1670/1671; YY/T0700 and so on.

### Technical date:

1. Measuring pressure points: 0kPa; 1.75kPa; 3.5kPa; 7kPa; 14kPa; 20kPa;
2. Pressing method: manual adjustment;
3. Timing range: 0 ~ 99.99min;
4. Sample size: 75mm × 75mm;
5. Pressure area: 28.27cm<sup>2</sup>;
6. Blocking net: 304 stainless steel; open space ≥50%; bending ≤5mm at 14kPa;
7. Clamp torque: 13.5N · m;
8. Air source: ≥0.4kg / cm<sup>2</sup>;
9. Power source: AC220V ± 10% 50HZ;
10. Power: 250W;
11. Net weight: 30Kg;
12. Dimensions: 450 × 400 × 550mm.



Tel: [+86-769-3365 7986](tel:+86-769-33657986)  
Website: [www.zltester.com](http://www.zltester.com)

Fax: [+86-769-3365 7987](tel:+86-769-33657987)  
E-mail: [zhongli@zltester.com](mailto:zhongli@zltester.com)